

IKA® Mixing
Magnetic stirrers with heating

hot



HS 7

	C-MAG HS 7
Technical data	
Stirring quantity (H ₂ O)	10 l
Motor rating input / output	15 / 1,5 W
Speed display	scale (0 - 6)
Speed range	100 – 1.500 min ⁻¹
Max. magnetic bar (L x Ø)	50 x 8mm
Heating function	
Heat output	1000 W
Heating plate temperature range	50 °C – 500 °CC
Heating plate setting accuracy	± 10 K
Safety circuit	(fixed) 550 °C
Control accuracy with sensor	ETS-D5 ± 0,5 K ETS-D6 ± 0,2 K
Heating plate	
Material	glass - ceramics
Dimensions	200 x 200
General data	
Dimensions (W x D x H)	220 x 330 x 105 mm
Weight	5 kg
Protection class acc. to DIN EN 60529	IP 21
Ident. No. (230 V / 115 V)	3581200 / 3581226