

BacT/ALERT 3D™
60



*Specific colorimetric
technology with visual
positive signal.*

Simply Smaller and Smarter

BacT/ALERT® 3D 60 is the latest addition to the BacT/ALERT 3D product range and is capable of blood culture, body fluid or mycobacterial testing. It confirms our commitment to providing continuous innovation to all laboratories.

Our space-saving design has been optimized still further. Combined with all the advantages of the larger volume modules, the BacT/ALERT® 3D 60 ensures all laboratories can benefit from this technology.


BacT/ALERT 3D 60 easily handles an annual workload of up to 3,600 blood/body fluid culture bottles, or up to 500 mycobacteria tests.

Providing you with a high level of convenience, BacT/ALERT 3D 60 simplifies your daily life.



Bringing the advantages of BacT/ALERT 3D to all laboratories

Unloaded Weight : 41Kg (90lbs)
Loaded Weight : 45Kg (98lbs)



Wand

- **Reduced risk of errors** with user-directed bottle loading and instant bottle recognition.
- **Immediate notification of results** due to the continuous bottle monitoring system in each cell.
- **Automatic quality control and instant bottle recognition** with the patented cell flag in each cell.
- **Peace-of-mind** through instant visual and audible alerts for positive bottles, as well as remote alarm capability.
- **Comprehensive data management** with an easy-to-use text-free interface.
- **Time, labor, and money saving** with one of the industry's lowest false-positive rates.



Load

Specifications

Capacity

- 60 cells

Electrical Power

- **Services Requirements**
- 100/240 Volts (50-60 HZ)
- 220/240 Volts (50-60 HZ)

Power Consumed in Watts

- 100/240 V.A.C. 72 Watts typical
- 220/240 V.A.C. 72 Watts typical

Heat Dissipated

- 245 Btu/Hr. 2840 maximum

Sound Emission

- 50 dB



Check

Contact your local bioMérieux representative to find out more about BacT/ALERT 3D 60.

bioMérieux S.A.
69280 Marcy l'Etoile
France
Tel. : (33) 04 78 87 20 00
Fax. : (33) 04 78 87 20 90

www.biomerieux.com

Manufactured by :
bioMérieux Inc.
100 Rodolphe Street,
Durham, NC 27712
USA
Tel. : 00 1 919 620 2000
Fax. : 00 1 919 620 2211

bioMérieux Inc. is a wholly-owned subsidiary of bioMérieux sa.

